

## Stirrer

### Hotplate Magnetic Stirrer

#### Features of BS-2H/3H:

- \* Cast Aluminum with surface spraying plate.
- \* Max. Temp.: 380°C, max. speed: 2000rpm.
- \* Heating and stirring can proceed simultaneously.
- \* Heating power and stirring speed is steplessly adjusted.



BS-2H

BS-3H

#### Features of BS-4H/4HC:

- \* Ceramic plate.
- \* PT100 sensor is standard for BS-4HC.
- \* With the adjustable rack and stirrer bar.
- \* Max. Temp.: 380°C, max. speed: 2000rpm.
- \* Heating and stirring can proceed simultaneously.
- \* Heating power and stirring speed is steplessly adjusted.
- \* Enclosed hot plate with flame protection, fast warm-up and durability.



BS-4H

BS-4HC

#### Technical Parameters:

Model	BS-2H	BS-3H	BS-4H	BS-4HC
Work Plate Size	120*120mm	170*170mm	190*190mm	
Plate Load(kg)	3	10		
Work Plate Material	Cast Aluminum with surface spraying		Ceramic	
Max. Stirring Capacity	2L	5L		
Speed Range	100~2000rpm		200~2000rpm	
Heating Temp	380°C			350 °C
Power Consumption	180W	500W	600W	
Power Supply	AC110/220V±10%, 50/60Hz			
Standard Accessory	Adjustable rack, stirrer bar; PT100 sensor (only for BS-4HC)			
External Size(W*D*H) mm	200*120*90	265*185*190	320*190*125	
Package Size(W*D*H) mm	250*180*120	280*210*220	390*240*230	
Net Weight(kg)	1.9	2.2	5	
Gross Weight(kg)	3	3.4	6	